

INHOUSE TESTING

Test Name	Test Code	Dept.	Acceptable Blood Tube(s)							Sample Type	Min. Volume	Storage/Stability			Special Instructions
			Light GRN LiHep	Dark GRN LiHep	Gold Silica	Red Silica	Blue NaCit	Lavender EDTA	Pink EDTA			Ambient (18-25°C)	Fridge (2-8°C)	Frozen (<-15°C)	
Acetaminophen Level	ACET	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Activated Partial Thromboplastin	PTT	HEME	NO	NO	NO	NO	[OK]	NO	NO	Plasma	0.5 mL	[4 hours]	N/A	14 days	Fill tube completely
Alanine Aminotransferase	ALT	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Albumin	ALB	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Alcohol (ETOH)	ETOH	CHEM	[OK]	NO	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	DO NOT STORE – RUN ASAP!			Fill tube to avoid aliquoting
Alkaline Phosphatase	ALKP	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Ammonia	NH3	CHEM	[OK]	NO	NO	NO	NO	NO	NO	Plasma	0.3 mL	30 min	30 min	N/A	Draw on ice, analyze w/o delay
Amylase	AMY	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Antenatal Rhogam	ANRR	BB	NO	NO	NO	NO	NO	OK	[OK]	Whole Blood	1 mL	N/A	[14 days]	N/A	No blood band
Antibody Screen	ABS	BB	NO	NO	NO	OK	NO	OK	[OK]	Plasma/Serum	0.5 mL	N/A	[14 days]	N/A	No blood band
Aspartate Aminotransferase	AST	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
BHCG (qualitative)	HCG	HEME	NO	NO	[OK]	OK	NO	NO	NO	Serum	0.1 mL	2 hours	[48 hours]	N/A	Female only
BHCG (quantitative)	BHCG	CHEM	OK	OK	[OK]	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	6 months	Female only
Bilirubin, Direct	DBIL	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	Protect from light
Bilirubin, Total	TBIL	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	Protect from light
Bilirubin, Total and Direct	TD	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	Protect from light
Blood Type	BT	BB	NO	NO	NO	NO	NO	[OK]	OK	Whole Blood	1 mL	N/A	[14 days]	N/A	No blood band
Blood Type, Confirmation	CONFBT	BB	NO	NO	NO	NO	NO	[OK]	OK	Whole Blood	1 mL	N/A	[14 days]	N/A	No blood band
Blood Urea Nitrogen	BUN	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
B-Type Natriuretic Peptide	BNP	CHEM	NO	NO	NO	NO	NO	[OK]	NO	Plasma	0.5 mL	7 hours	[7 hours]	7 days	
Calcium, Ionized	ICA	CHEM	NO	[OK]	NO	NO	NO	NO	NO	Whole Blood	0.5 mL	[1 hour]	N/A	N/A	
Calcium, Total	CA	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Carbon Dioxide	CO2	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Carboxyhemoglobin Panel	COX	CHEM	NO	[OK]	NO	NO	NO	NO	NO	Whole Blood	0.5 mL	14 days	[14 days]	N/A	
Carcinoembryonic Antigen	CEA	CHEM	OK	OK	[OK]	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	7 days	
Chem-8 Profile (BMP)	CX8P	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Chloride	CL	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Cholesterol, Total	CHOL	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Coagulation Survey	CS	HEME	NO	NO	NO	NO	REQ	REQ	NO	Plasma Whole Blood	1.5 mL 0.3 mL	[4 hours]	N/A	N/A	Blue + Lavender top required
Complete Blood Count	CBC	HEME	NO	NO	NO	NO	NO	[OK]	OK	Whole Blood	0.3 mL	[24 hours]	48 hours	N/A	

Comprehensive Metabolic Panel	CMP	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Cord Blood Evaluation	CORDB	BB	NO	NO	NO	NO	NO	[OK]	NO	Cord Blood	1 mL	N/A	[14 days]	N/A	Floor collected specimen
Cortisol	CORT CORST	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	7 days	
C-Reactive Protein	CRP	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Creatine Kinase	CK	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	4 hours	[12 hours]	48 hours	
Creatinine	CRET	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
D-Dimer	DIM	HEME	NO	NO	NO	NO	[OK]	NO	NO	Plasma	0.5 mL	[8 hours]	14 days	6 months	Fill tube completely
DIC Screen	DIC	HEME	NO	NO	NO	NO	REQ	REQ	NO	Plasma Whole Blood	1.5 mL 1 mL	[4 hours]	N/A	N/A	Blue + Lavender top required
Digoxin Level	DIG	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Direct Antiglobulin Test	DAT	BB	NO	NO	NO	NO	NO	[OK]	OK	Whole Blood	1 mL	72 hours	[72 hours]	N/A	No blood band
Electrolyte Profile	LYTE	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Estradiol, Sensitive	ESTRDL	CHEM	[OK]	NO	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[7 days]	5 days	
Ferritin	FERR	CHEM	[OK]	NO	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	7 days	
Fibrinogen	FIBRIN	HEME	NO	NO	NO	NO	[OK]	NO	NO	Plasma	0.5 mL	[8 hours]	N/A	14 days	Fill tube completely
Folate	FOL	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	N/A	[48 hours]	6 months	
Follicle Stimulating Hormone	FSH	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	7 days	
Gamma-Glutamyl Transferase	GGT	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Gentamicin Level	GENR GENT GENP	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	Time specific
Glucose	GLU	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Glycohemoglobin	GLYC	CHEM	NO	NO	NO	NO	NO	[OK]	NO	Whole Blood	0.3 mL	6 days	[30 days]	45 days	
HDL Cholesterol	1HDL	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/serum	0.3 mL	8 hours	[48 hours]	1 month	
Hemoglobin & Hematocrit	HNHL	HEME	NO	NO	NO	NO	NO	[OK]	OK	Whole Blood	0.3 mL	[24 hours]	48 hours	N/A	
Hemogram	HGML	HEME	NO	NO	NO	NO	NO	[OK]	OK	Whole Blood	0.3 mL	[24 hours]	48 hours	N/A	
Heparin, Low-Molecular Weight	LMWH	HEME	NO	NO	NO	NO	[OK]	NO	NO	Plasma	0.5 mL	[2 hours]	N/A	N/A	Fill tube completely
Heparin, Unfractionated	1UFH	HEME	NO	NO	NO	NO	[OK]	NO	NO	Plasma	0.5 mL	[2 hours]	N/A	N/A	Fill tube completely
Hepatic Function Panel	HFP	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma	0.5 mL	8 hours	[48 hours]	1 month	
Hepatitis A IgM Antibody	HAVAB	CHEM	OK	NO	[OK]	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[7 days]	1 month	
Hepatitis B Core IgM Antibody	HBCAB	CHEM	OK	NO	[OK]	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[7 days]	1 month	
Hepatitis B Surface Antibody	HBSAB	CHEM	NO	NO	[OK]	OK	NO	NO	NO	Serum	0.5 mL	8 hours	[48 hours]	1 month	
Hepatitis B Surface Antigen	HBSAG	CHEM	OK	NO	[OK]	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Hepatitis C IgG Antibody	HEPC	CHEM	OK	NO	[OK]	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	For ages >10
Hepatitis Profile	HEPP	CHEM	OK	NO	[OK]	OK	NO	NO	NO	Plasma/Serum	1 mL	8 hours	[48 hours]	1 month	For ages >10
HIV 1+2 Antibody	HIVCOM	CHEM	OK	NO	[OK]	NO	NO	REFLEX	NO	Serum		8 hours	[48 hours]	1 month	For ages >2 Reflex needs lavender top
Infectious Mononucleosis	MONO	BB	NO	NO	[OK]	OK	NO	OK	NO	Serum	0.5 mL	N/A	[24 hours]	N/A	

Iron, Total	FET	CHEM	OK	NO	[OK]	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Iron/Total Iron Binding Capacity	FE	CHEM	OK	NO	[OK]	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Kleihauer-Betke	KLBT	BB	NO	NO	NO	NO	NO	OK	[OK]	Whole Blood	1 mL	N/A	[14 days]	N/A	
Lactate Dehydrogenase	LD	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	[8 hours]	N/A	N/A	
Lactic Acid	LACT	CHEM	NO	[OK]	NO	NO	NO	NO	NO	Whole Blood	0.5 mL	30 min	N/A	N/A	Analyze without delay
LDL Cholesterol, Direct	DLDL	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Lipase	LIPA	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	4 hours	[48 hours]	1 month	
Lipid Profile	LIPID	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	12 hours fasting
Lithium Level	LI	CHEM	NO	NO	[OK]	OK	NO	NO	NO	Serum	0.5 mL	8 hours	[7 days]	1 month	
Luteinizing Hormone	LH	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	6 months	
Magnesium	MG	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Osmolality	OSM	CHEM	NO	NO	[OK]	OK	NO	NO	NO	Serum	0.5 mL	24 hours	[48 hours]	7 days	
Parathyroid Hormone, Intact	PTHIO	CHEM	[OK]	OK	OK	OK	NO	OK	NO	Plasma/Serum	0.5 mL	4-8 hours	[8-48 hours]	6 months	Serum has a shorter stability time
Parathyroid Hormone, Intact w/ Calcium	PTHIM	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	4-8 hours	[8-48 hours]	6 months	Serum has a shorter stability time
pH, Venous	VPH	CHEM	NO	[OK]	NO	NO	NO	NO	NO	Whole Blood	0.5 mL	30 min	N/A	N/A	Analyze without delay
Phenobarbital Level	PBARB	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Phenytoin (Dilantin) Level	DILN	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Phosphorus	PHOS	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Platelet Function Analysis	PFA	HEME	NO	NO	NO	NO	[OK] x2	NO	NO	Whole Blood	5.4 mL	[4 hours]	N/A	N/A	Do not send through tube system
Post Delivery Rhogam	PDR	BB	NO	NO	NO	NO	NO	[OK]	OK	Whole Blood	0.5 mL	5 days	[5 days]	N/A	
Potassium	K	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Prealbumin	PALB	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Procalcitonin	PROCL	CHEM	NO	NO	[OK]	OK	NO	NO	NO	Serum	0.5 mL	16 hours	[48 hours]	75 days	
Progesterone	PROG	CHEM	NO	NO	NO	[OK]	NO	NO	NO	Serum	0.5 mL	8 hours	[48 hours]	7 days	
Prolactin	PROL	CHEM	OK	NO	[OK]	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	7 days	
Protein, Total	TP	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Protim/INR	PTINR	HEME	NO	NO	NO	NO	[OK]	NO	NO	Plasma	0.5 mL	[8 hours]	N/A	14 days	Fill tube completely
Rapid Plasma Reagin Titer	RPR	BB	NO	NO	[OK]	OK	NO	OK	NO	Plasma/Serum	0.2 mL	Not Established			
Renal Profile	RENP	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Reticulocyte Count	RETLH	HEME	NO	NO	NO	NO	NO	[OK]	OK	Whole Blood	0.5 mL	[24 hours]	48 hours	N/A	
Rubella IgG Antibody	RUBE	CHEM	NO	NO	[OK]	OK	NO	NO	NO	Serum	0.5 mL	8 hours	[48 hours]	1 year	
Salicylate Level	SALI	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Sodium	NA	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
T3, Free	FT3	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	7 days	
T4, Free	FT4	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	7 days	

Tegretol (Carbamazepine) Level	TEG	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Theophylline Level	THEO	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Thyroid Stimulating Hormone	TSH TSHR	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	18 hours	[7 days]	3 months	
Tobramycin Level	TOBR TOBP TOBT	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	Time specific
Transferrin	TRANSF	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Triglycerides	TG	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	
Troponin-I	TROP	CHEM	[OK]	OK	NO	NO	NO	NO	NO	Plasma	0.5 mL	2 hours	[24 hours]	6 months	
Type and Screen	TYSC	BB	NO	NO	NO	OK	NO	OK	REQ	Whole Blood	1 mL	N/A	[14 days]	N/A	Blood band required
Type and Screen, Preadmit	PRTYSC	BB	NO	NO	NO	NO	NO	OK	REQ	Whole Blood	1 mL	N/A	[14 days]	N/A	Blood band required
Type and Screen, Prenatal	TYSCP	BB	NO	NO	NO	NO	NO	OK	[OK]	Whole Blood	1 mL	N/A	[14 days]	N/A	No blood band
Uric Acid Level	URIC	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Valproic Acid (Depakote) Level	VPA	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.3 mL	8 hours	[48 hours]	1 month	
Vancomycin Level	VANR VANP VANT	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[48 hours]	1 month	Time specific
Vitamin B12 Level	B12	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	8 hours	[24 hours]	6 months	
Vitamin D (1,25-OH) Level	VITD	CHEM	[OK]	OK	OK	OK	NO	NO	NO	Plasma/Serum	0.5 mL	72 hours	[7 days]	1 year	
Westergren Sed Rate	WSR	HEME	NO	NO	NO	NO	NO	[OK]	OK	Plasma/Serum	0.5 mL	[12 hours]	24 hours	N/A	

[...] = Preferred tube type/storage | REQ = Required tube |

SENDOUT TESTING

Test Name	Test Code	Acceptable Blood Tube(s)								Sample Type	Min. Volume	Storage/Stability			Special Instructions
		Light GRN LiHep	Dark GRN LiHep	Gold Silica	Red Silica	Blue NaCit	Lavender EDTA	Pink EDTA	Other			Ambient (18-25°C)	Fridge (2-8°C)	Frozen (<-15°C)	
1,25-Dihydroxyvitamin D	DHVD	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	1.5 mL	7 days	[30 days]	30 days	Aliquot to new container
17-Hydroxyprogesterone	OHPG	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.25 mL	7 days	[28 days]	28 days	Aliquot to new container
Acetylcholine Receptor (Muscle AChR) Binding Antibody	FABAB	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	1 mL	72 hours	[28 days]	28 days	Aliquot to new container
Adrenocorticotrophic Hormone	ACTHM	NO	NO	NO	NO	NO	[OK]	OK	NO	Plasma	0.75 mL	2 hours	3 hours	[28 days]	Refer to test catalog!
Aldolase	ALDOM	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.5 mL	N/A	[7 days]	60 days	Centrifuge/aliquot within 1 hour
Aldosterone	ALDS	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	1.2 mL	4 days	[28 days]	30 days	Refer to test catalog!
Alpha-1-Antitrypsin	AAT	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.5 mL	28 days	[28 days]	28 days	Aliquot to new container
Alpha-Fetoprotein (AFP), Single Marker Screen, Maternal	MAFP1	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	7 days	[14 days]	90 days	Refer to test catalog!
Alpha-Fetoprotein (AFP) Tumor Marker	AFPTM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	N/A	[7 days]	90 days	Aliquot to new container
Amino Acids, Quantitative	AAQP	NO	NO	NO	NO	NO	NO	NO	[Na Hep]	Plasma	0.3 mL	N/A	N/A	[14 days]	Refer to test catalog!
Angiotensin Converting Enzyme	ACEM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	24 hours	[7 days]	180 days	Aliquot to new container
Antibody to Extractable Nuclear Antigen Evaluation	ENAE	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.35 mL	N/A	[21 days]	21 days	Aliquot to new container
Antimullerian Hormone	AMH1	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	7 days	[7 days]	180 days	Aliquot to new container
Antineutrophil Cytoplasmic Antibodies Vasculitis Panel	VASC	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	N/A	[21 days]	21 days	Aliquot to new container
Antinuclear Antibodies (ANA)	ANA2M	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Antistrep-O Titer	ASOM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	7 days	[28 days]	28 days	Aliquot to new container
Antithrombin Activity	ATTE	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Antithrombin Antigen	ATTI	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Ascorbic Acid (Vitamin C)	VITCM	OK	[OK]	NO	NO	NO	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Collect on ice and protect from light Aliquot to amber tube and freeze
Aspergillus (Galactomannan) Antigen	ASPAGS	NO	NO	[OK]	NO	NO	NO	NO	NO	Serum	1 mL	N/A	[14 days]	14 days	Avoid exposure to air Send in original tube
BCR/ABL1 Qualitative w/ Reflex to BCR/ABL1 p190 or p210 Quantitative	BCRFLX	NO	NO	NO	NO	NO	OK	[OK]	ACD	Whole Blood	8 mL	72 hours	[72 hours]	N/A	Send in original tube Label specimen as blood

BCR/ABL1, p190, mRNA Detection, RT-PCR, Quantitative, Monitoring Assay	BA190	NO	NO	NO	NO	NO	OK	[OK]	ACD	Whole Blood	8 mL	72 hours	[5 days]	N/A	Send in original tube Label specimen as blood
BCR/ABL1, p210, mRNA Detection, RT-PCR, Quantitative, Monitoring CML	BA210	NO	NO	NO	NO	NO	OK	[OK]	ACD	Whole Blood	4 mL	72 hours	[72 hours]	N/A	Send in original tube Label specimen as blood
Beta-2 Glycoprotein 1 Antibodies, IgG and IgM	B2GMG	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Beta-2-Microglobulin	B2M	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	14 days	[28 days]	28 days	Aliquot to new container
Beta-Human Chorionic Gonadotropin, Quantitative	BHCSTM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	7 days	[7 days]	90 days	Aliquot to new container
Beta-Hydroxybutyrate	BHYD	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.25 mL	N/A	7 days	[7 days]	Centrifuge/aliquot within 2 hours
Bile Acids, Fractionated and Total	BAFS	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.3 mL	90 days	[90 days]	90 days	12-14 hours fasting required
Bile Acids, Total	BILEA	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.25 mL	24 hours	[7 days]	30 days	Refer to test catalog!
Breast Carcinoma-Associated Antigen	C2729M	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.25 mL	4 days	[7 days]	90 days	Aliquot within 2 hours
Calcitonin	CATN	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	8 hours	24 hours	[90 days]	Collect on ice Use chilled centrifuge cup
Cancer Antigen 125 (CA 125)	CA25	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	N/A	[5 days]	168 days	Aliquot to new container
Cancer Antigen 15-3 (CA 15-3)	CA153M	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	N/A	[7 days]	90 days	Aliquot to new container
Candida albicans (Monilia), IgE	CDAB	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	N/A	[14 days]	90 days	Aliquot to new container
Carbohydrate Antigen 19-9 (CA 19-9)	CA19	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	8 hours	[7 days]	90 days	Centrifuge within 2 hours
Celiac Disease Serology Cascade	CDSP	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	2 mL	N/A	[14 days]	21 days	Aliquot to new container
Centromere Antibodies, IgG	CMAM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.35 mL	N/A	[21 days]	21 days	Aliquot to new container
Ceruloplasmin	CERUM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	N/A	[7 days]	30 days	Aliquot to new container
Chromogranin A	CGAK	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.2 mL	48 hours	48 hours	[90 days]	Aliquot to new container
Chromosome Analysis, Congenital Disorders	CHRCB	NO	NO	NO	NO	NO	NO	NO	[Na Hep]	Whole Blood	2 mL	Not Specified			Refer to test catalog!
Chromosome Microarray, Congenital	CMACB	NO	NO	NO	NO	NO	[REQ]	NO	[REQ]	Whole Blood	2 mL	Not Specified			Refer to test catalog!
Clozapine	CLZ	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.6 mL	7 days	[28 days]	28 days	Refer to test catalog!
Coagulation Factor II Assay	FACT2	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Coagulation Factor V Assay	FACTVR	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Coagulation Factor VII Assay	FACT7	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Coagulation Factor VIII Activity Assay	F8A	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Coagulation Factor IX Assay	FACT9	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Coagulation Factor X Assay	FACTX	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Coagulation Factor XI Assay	FACT11	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!

Coagulation Factor XII Assay	FACT12	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Coccidioides Antibody, Complement Fixation and Immunodiffusion	SCOC	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	1.2 mL	N/A	[14 days]	14 days	Aliquot to new container
Complement C3	C3M	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	72 hours	[7 days]	28 days	Aliquot to new container
Complement C4	C4M	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	72 hours	[7 days]	28 days	Aliquot to new container
Complement, Total	COMPT	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	N/A	N/A	[28 days]	Collect on ice Immediately freeze after centrifuging
Connective Tissue Disease Cascade	CTDC	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.7 mL	N/A	[21 days]	21 days	Aliquot to new container
Copper	CUS	NO	NO	NO	NO	NO	NO	NO	[Royal Blue]	Serum	0.2 mL	28 days	[28 days]	28 days	Refer to test catalog!
Coxiella burnetii (Q Fever), Molecular Detection, PCR	CBSRP	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.5 mL	N/A	[7 days]	7 days	Centrifuge/aliquot within 2 hours
C-Peptide	CPR	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	N/A	7 days	[30 days]	Refer to test catalog!
C-Reactive Protein, High Sensitivity	HSCRPM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.2 mL	N/A	[7 days]	30 days	Aliquot to new container
Cryptococcus Antigen Screen with Titer	SLFA	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	N/A	[14 days]	14 days	Aliquot to new container
Cryptococcus Antigen Titer, Lateral Flow Assay	SLFAT	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.3 mL	N/A	[14 days]	14 days	Aliquot to new container
Cyclic Citrullinated Peptide Antibodies, IgG	CCPM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Cyclosporine	CYSRP	NO	NO	NO	NO	NO	[OK]	NO	NO	Whole Blood	1 mL	14 days	[14 days]	14 days	Send in original tube
Cytomegalovirus (CMV) Antibodies, IgM and IgG	CMVAB	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.8 mL	N/A	[14 days]	14 days	Aliquot to new container
Cytomegalovirus (CMV) DNA Detection and quantification by RT-PCR	CMVQU	NO	NO	NO	NO	NO	[OK]	OK	NO	Plasma	0.6 mL	N/A	6 days	[84 days]	Centrifuge/aliquot within 2 hours
Cytoplasmic Neutrophil Antibodies	MANCA	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Dehydroepiandrosterone (DHEA)	DHEAM	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.5 mL	6 hours	21 days	[28 days]	Aliquot to new container
Dehydroepiandrosterone Sulfate	DHESM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	7 days	[14 days]	30 days	Aliquot to new container
Dilute Russell Viper Venom Time (DRVVT), with reflex	DRVTI	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
DNA Double-Stranded (dsDNA) Antibodies by <i>Critidia luciliae</i> , IgG	CRITH	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.3 mL	N/A	[21 days]	21 days	Aliquot to new container
DNA Double-Stranded (dsDNA) Antibodies with reflex, IgG	ADNAR	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.3 mL	N/A	[21 days]	21 days	Aliquot to new container
Endomysial Antibodies, IgA	EMAM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	14 days	[14 days]	30 days	Aliquot to new container
Enterovirus, Molecular Detection, PCR	ENTP	NO	NO	NO	NO	NO	[OK]	OK	NO	Plasma	0.3 mL	N/A	[7 days]	7 days	Aliquot to new container
Epstein-Barr Virus Antibody Profile	EBVAB	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	48 hours	[14 days]	Aliquot to new container
Epstein-Barr Virus DNA Detection and Quantification	EBVQU	NO	NO	NO	NO	NO	[OK]	OK	NO	Plasma	0.5 mL	N/A	6 days	[84 days]	Centrifuge/aliquot within 2 hours
Erythropoietin	EPO	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	7 days	[14 days]	14 days	Morning collection preferred
Estrone	ESTRM	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.7 mL	28 days	[28 days]	28 days	Centrifuge/aliquot within 2 hours

Everolimus	EVROL	NO	NO	NO	NO	NO	[OK]	OK	NO	Whole Blood	1 mL	14 days	[28 days]	28 days	Send in original tube
Factor V Leiden (R506Q) Mutation	FACVL	NO	NO	NO	NO	OK	[OK]	OK	ACD (B)	Whole Blood	1 mL	[14 days]	14 days	14 days	Send in original tube
Gastrin	GSTRN	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	N/A	24 hours	[30 days]	Refer to test catalog!
Gliadin (Deamidated) Antibodies Evaluation, IgG and IgA	DGLDN	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Gliadin (Deamidated) Antibody, IgA	DAGL	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Gliadin (Deamidated) Antibody, IgG	DGGL	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Glomerular Basement Membrane Antibodies, IgG	GBMAB	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.35 mL	N/A	[21 days]	21 days	Aliquot to new container
Glucose 6-Phosphate Dehydrogenase	G6PDM	NO	NO	NO	NO	NO	OK	OK	[ACD]	Whole Blood	1 mL	N/A	[20 days]	N/A	Send in original tube
Glutamic Acid Decarboxylase (GAD65)	GAD65	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	1mL	72 hours	[28 days]	28 days	Aliquot to new container
Hantavirus Antibody (IgG, IgM)	FHVGM	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.5 mL	7 days	[14 days]	30 days	Aliquot to new container
Haptoglobin	HAPTM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	14 days	[28 days]	28 days	Aliquot to new container
Hepatitis A IgM Antibody	HAIGM	NO	NO	[OK]	NO	NO	NO	NO	NO	Serum	0.6 mL	72 hours	6 days	[90 days]	Centrifuge/aliquot within 2 hours
Hepatitis B Core Total Antibodies	HBCM	NO	NO	[OK]	NO	NO	NO	NO	NO	Serum	0.6 mL	72 hours	6 days	[90 days]	Centrifuge/aliquot within 2 hours
Hepatitis B Surface Antigen	HBAG	NO	NO	[OK]	NO	NO	NO	NO	NO	Serum	0.7 mL	72 hours	6 days	[90 days]	Centrifuge/aliquot within 2 hours
Hepatitis B Virus (HBV) DNA Detection and Quantification by RT-PCR	HBVQU	NO	NO	[OK]	NO	NO	NO	NO	NO	Serum	0.8 mL	N/A	6 days	[84 days]	Centrifuge/aliquot within 2 hours
Hepatitis C RNA Detection and Quantification by RT-PCR	HCVQN	NO	NO	[OK]	NO	NO	NO	NO	NO	Serum	0.8 mL	N/A	6 days	[84 days]	Centrifuge/aliquot within 2 hours
Hepatitis C RNA Quantification w/ Reflex to HCV Genotype	HCVQG	NO	NO	[OK]	NO	NO	NO	NO	NO	Serum	1.6 mL	N/A	72 hours	[42 days]	Centrifuge/aliquot within 2 hours
Hepatitis C Virus Antibody w/ Reflex to HCV RNA, PCR, Symptomatic	HCVDX	NO	NO	[OK]	NO	NO	NO	NO	NO	Serum	0.9 mL	N/A	6 days	[84 days]	Centrifuge/aliquot within 2 hours
Hepatitis C Virus Genotype	HCVG	NO	NO	[OK]	NO	NO	NO	NO	NO	Serum	1.5 mL	N/A	72 hours	[42 days]	Centrifuge/aliquot within 2 hours
Hereditary Hemochromatosis, HFE Variant Analysis	HFE	OK	OK	NO	NO	OK	[OK]	OK	ACD	Whole Blood	0.5 mL	[4 days]	14 days	N/A	Send in original tube
HIV-1/HIV-2 Antibody Confirmation and Differentiation	HIVDI	NO	NO	NO	NO	NO	[OK]	OK	NO	Plasma	0.8 mL	N/A	6 days	[30 days]	Centrifuge/aliquot within 2 hours
HIV-1/HIV-2 RNA Detection	HIP12	NO	NO	NO	NO	NO	[OK]	OK	NO	Plasma	0.8 mL	N/A	6 days	[40 days]	Centrifuge/aliquot within 2 hours
HIV-1 RNA Detection and Quantification	HIVQN	NO	NO	NO	NO	NO	[OK]	OK	NO	Plasma	0.8 mL	N/A	6 days	[84 days]	Centrifuge/aliquot within 2 hours
HLA-B27	LY27B	NO	NO	NO	NO	NO	OK	[OK]	NO	Whole Blood	1 mL	[4 days]	N/A	N/A	Send in original tube
Homocysteine, Total	HCYSP	OK	OK	NO	NO	NO	[OK]	NO	NO	Plasma	0.3 mL	28 days	[28 days]	309 days	Draw on ice Centrifuge/aliquot within 4 hours
HSV Type 1- and Type 2-Specific Antibodies, IgG	MHSVG	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[14 days]	14 days	Aliquot to new container
HSV, Molecular Detection, PCR	LHSV B	NO	NO	NO	NO	NO	[OK]	NO	NO	Whole Blood	0.4 mL	N/A	[7 days]	7 days	Send in original tube

Immunofixation Heavy Chain Type Delta and Epsilon	IFXD	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	1 mL	14 days	[28 days]	28 days	Aliquot to new container
Immunoglobulin A (IgA)	MIGA	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	14 days	[28 days]	28 days	Aliquot to new container
Immunoglobulin D (IgD)	MIGD	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	21 days	[28 days]	28 days	Aliquot to new container
Immunoglobulin E (IgE)	MIGE	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	N/A	[14 days]	90 days	Aliquot to new container
Immunoglobulin Free Light Chains	FLCP	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	72 hours	[28 days]	28 days	Aliquot to new container
Immunoglobulin G (IgG)	MIGG	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	14 days	[28 days]	28 days	Aliquot to new container
Immunoglobulin M (IgM)	MIGM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	14 days	[28 days]	28 days	Aliquot to new container
Immunoglobulins (IgG, IgA, and IgM)	IMMG	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	14 days	[28 days]	28 days	Aliquot to new container
Infliximab Quantitation w/ Reflex to Antibodies to Infliximab	INFXR	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.5 mL	N/A	[28 days]	28 days	Refer to test catalog!
Inhibin A, Tumor Marker	INHA	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	7 days	[14 days]	90 days	Aliquot to new container
Inhibin B	INHB	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	7 days	[14 days]	90 days	Aliquot to new container
Insulin	INSM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	N/A	7 days	[180 days]	Refer to test catalog!
Insulin Antibodies	INAB	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	1 mL	72 hours	[28 days]	28 days	Aliquot to new container
Insulin-Like Growth Factor-1, Mass Spectrometry	IGFMS	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.3 mL	72 hours	72 hours	[28 days]	Aliquot to new container
Intrinsic Factor Blocking Antibody	IFBA	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	N/A	[14 days]	14 days	Refer to test catalog!
JAK2 V617F Mutation Detection	JAK2B	NO	NO	NO	NO	NO	[OK]	OK	ACD	Whole Blood	1 mL	[7 days]	7 days	N/A	Send in original tube
Lamotrigine	LAMO	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.5 mL	28 days	[28 days]	28 days	Centrifuge/aliquot immediately
Lead, Capillary, w/ Demographics	PBDC	NO	NO	NO	NO	NO	NO	NO	[Micro MAP]	Whole Blood	0.1 mL	28 days	[28 days]	28 days	Refer to test catalog!
Lead, Venous, w/ Demographics	PBDV	NO	NO	NO	NO	NO	NO	NO	[Royal Blue]	Whole Blood	0.1 mL	28 days	[28 days]	28 days	Refer to test catalog!
Levetiracetam	LEVE	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.5 mL	28 days	[28 days]	28 days	Centrifuge/aliquot within 2 hours
Liver/Kidney Microsome Type 1 Antibodies	LKM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Lupus Anticoagulant Profile	ALUPP	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	3 mL*	N/A	N/A	[14 days]	Refer to test catalog!
Lyme Disease Antibody, Immunoblot	LYWB	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	N/A	[14 days]	30 days	Aliquot to new container
Lyme Disease Serology	LYMEM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[10 days]	30 days	Aliquot to new container
Lysozyme (Muramidase)	MURA	NO	NO	NO	NO	NO	[OK]	OK	NO	Plasma	1 mL	N/A	N/A	[21 days]	Centrifuge/aliquot within 2 hours
Measles (Rubeola) Antibodies, IgG	ROPG	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[14 days]	14 days	Aliquot to new container
Methotrexate	MTHX	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.25 mL	72 hours	[14 days]	28 days	Protect from light Aliquot to amber tube within 2 hours
Methylmalonic Acid, Quantitative	MMAS	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	48 days	[48 days]	48 days	Aliquot to new container

Mitochondrial Antibodies (M2)	MITAB	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Monoclonal Protein Study, Quantitative	QMPSS	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	1.5 mL	7 days	[28 days]	28 days	Refer to test catalog!
MPN JAK2 V617F w/ Reflex to CALR and MPL	MPNR	NO	NO	NO	NO	NO	[OK]	OK	ACD	Whole Blood	0.5 mL	[7 days]	7 days	N/A	Send in original tube
Mumps Virus Antibody, IgG	MPPG	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[14 days]	14 days	Aliquot to new container
Mycophenolic Acid	MPNA	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.25 mL	21 days	[28 days]	28 days	Aliquot to new container
Newborn Aneuploidy Detection, FISH	NADF	NO	NO	NO	NO	NO	OK	OK	[Na Hep]	Whole Blood	1 mL	Not Specified			Refer to test catalog!
Nicotine and Metabolites	NICOS	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.5 mL	28 days	[28 days]	28 days	Centrifuge/aliquot within 2 hours
Parathyroid Hormone-Related Peptide	PTHPEP	NO	NO	NO	NO	NO	[OK]	OK	NO	Plasma	0.25 mL	N/A	N/A	[30 days]	Collect on ice Use chilled centrifuge cup
Parvovirus B19 Antibodies, IgG and IgM	PARVS	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[14 days]	14 days	Aliquot to new container
Phospholipid (Cardiolipin) Antibodies, IgA	ACLIPM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Phospholipid (Cardiolipin) Antibodies, IgG and IgM	CLPMG	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Protein C Activity	PACT	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Protein C Antigen	PCAG	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[2 years]	Refer to test catalog!
Protein S Activity	PSACT	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Protein S Antigen	PSTF	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Prothrombin G20210A Mutation	PTNT	NO	NO	NO	NO	OK	[OK]	OK	ACD (B)	Whole Blood	1 mL	[14 days]	14 days	14 days	Send in original tube
Quad Screen (Second Trimester) Maternal	QUAD1	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	7 days	[7 days]	90 days	Refer to test catalog!
QuantiFERON-TB Gold Plus	QFGM	NO	NO	NO	NO	NO	NO	NO	[TB Gold Kit]	Whole Blood	4 mL	N/A	[28 days]	N/A	Refer to test catalog!
Quantitative Lymphocyte Subsets: T, B and NK Cells	TBBS	NO	NO	NO	NO	NO	[OK]	OK	NO	Whole Blood	1 mL	[52 hours]	N/A	N/A	Send in original tube
Rapid Malaria/Babesia Smear	MLARIA	NO	NO	NO	NO	NO	[OK]	OK	NO	Whole Blood	0.5 mL	N/A	[Store & Transport]	N/A	***Requires 4 blood smears*** Refer to test catalog!
Renin Activity	PRAM	NO	NO	NO	NO	NO	[OK]	OK	NO	Plasma	1.15 mL	N/A	N/A	[14 days]	Refer to test catalog!
Reptilase Time	RPTL	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Respiratory Profile, Inland Northwest (OR, Central/Eastern WA)	INWAP	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	1.4 mL	N/A	[14 days]	90 days	Aliquot to new container
Rheumatoid Factor	RHUT	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	14 days	[14 days]	14 days	Aliquot to new container
Ribosome P Antibodies, IgG	RIB	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.35 mL	N/A	[21 days]	21 days	Aliquot to new container
Serotonin	SERWB	NO	NO	NO	NO	NO	[OK]	OK	NO	Whole Blood	1.5 mL	N/A	N/A	[90 days]	Refer to test catalog!
Sex Hormone-Binding Globulin	SHBGM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.5 mL	7 days	[14 days]	90 days	Aliquot to new container
Sirolimus	SIIRO	NO	NO	NO	NO	NO	[OK]	OK	NO	Whole Blood	1 mL	28 days	[28 days]	28 days	Send in original tube

Smooth Muscle Antibody Screen	SMA	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Syphilis IgG with Reflex, Enzyme Immunoassay	SYPHT	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	1 mL	N/A	N/A	[14 days]	Aliquot to new container
T3, Reverse	RT3M	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.4 mL	28 days	[28 days]	28 days	Aliquot to new container
T3, Total	T3M	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	N/A	[7 days]	30 days	Aliquot to new container
T4, Total Only	TT4M	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.75 mL	N/A	[7 days]	30 days	Aliquot to new container
Tacrolimus	TAKRO	NO	NO	NO	NO	NO	[OK]	OK	NO	Whole Blood	1 mL	14 days	[14 days]	14 days	Send in original tube
Testosterone, Total and Free	TGRP	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	1 mL	N/A	[14 days]	60 days	Aliquot to new container
Testosterone, Total, Mass Spectrometry	TTST	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.25 mL	N/A	[14 days]	60 days	Aliquot to new container
Thiamine (Vitamin B1)	VITB1M	NO	NO	NO	NO	NO	[OK]	OK	NO	Whole Blood	0.5 mL	N/A	N/A	[28 days]	Protect from light Aliquot to amber tube and freeze
Thrombin Time	TT	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Thyroglobulin Antibody	TGAB	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.5 mL	7 days	[7 days]	30 days	Centrifuge/aliquot within 2 hours
Thyroglobulin Mass Spectrometry	TGBMS	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.75 mL	72 hours	[7 days]	416 days	Aliquot to new container
Thyroglobulin, Tumor Marker Reflex	HTGR	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	1.25 mL	72 hours	[7 days]	30 days	Aliquot to new container
Thyroid Autoantibodies Profile	TAB	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.5 mL	7 days	[7 days]	30 days	Centrifuge/aliquot within 2 hours
Thyroperoxidase Antibodies	MTPO	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.5 mL	7 days	[7 days]	30 days	Centrifuge/aliquot within 2 hours
Tissue Transglutaminase Antibody, IgA	TTGA	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
Tissue Transglutaminase Antibody, IgG	TTGGM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[21 days]	21 days	Aliquot to new container
ToRCH Profile IgG	TCHG	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	1.2 mL	N/A	[14 days]	14 days	Aliquot to new container
Toxoplasma gondii Antibody, IgM and IgG	TXMGP	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.8 mL	N/A	[14 days]	14 days	Aliquot to new container
Tryptase	TRYPTM	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.2 mL	N/A	7 days	[14 days]	Centrifuge/aliquot immediately
Varicella-Zoster Antibody, IgG	VZPG	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[14 days]	14 days	Aliquot to new container
Viscosity	VISC	NO	NO	NO	[OK]	NO	NO	NO	NO	Serum	0.75 mL	14 days	[28 days]	28 days	Keep specimen at 37°C until centrifuged and aliquoted
Vitamin A	VTMNA	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.25 mL	28 days	[28 days]	28 days	Centrifuge/aliquot within 2 hours
Vitamin E	VITEM	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.25 mL	7 days	[44 days]	44 days	Protect from light Aliquot to amber tube within 24 hours
Vitamin K1	VITK1	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	0.75 mL	30 days	[30 days]	30 days	Centrifuge/aliquot within 2 hours
Von Willebrand Disease Profile	VWP	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	2 mL	N/A	N/A	[14 days]	Refer to test catalog!
Von Willebrand Factor Activity	VWACT	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!
Von Willebrand Factor Antigen	VWA	NO	NO	NO	NO	[OK]	NO	NO	NO	Plasma	0.5 mL	N/A	N/A	[14 days]	Refer to test catalog!

West Nile Virus Antibody, IgG and IgM	WNS	NO	NO	[OK]	OK	NO	NO	NO	NO	Serum	0.4 mL	N/A	[14 days]	14 days	Aliquot to new container
Zinc	ZNS	NO	NO	NO	NO	NO	NO	NO	[Royal Blue]	Serum	0.2 mL	28 days	[28 days]	28 days	Refer to test catalog!
Zinc Transporter 8 (ZnT8) Antibody	EZNT8	NO	NO	OK	[OK]	NO	NO	NO	NO	Serum	1 mL	72 hours	[28 days]	28 days	Aliquot to new container

[...] = Preferred tube type/storage | REQ = Required tube |